

Introduction

SimplePCI Automated Image Capture-Advanced Devices (AIC-AD) enhances **SimplePCI** to provide cost-effective control over all major motorized microscopes and select scanning confocal attachments. A Scan Wizard simplifies collection of image sequences for XYZ and wavelength over time (time-lapse). Image sequences are easily viewed in a data file for later analysis.

Expand the functionality of **SimplePCI** and **AIC-AD** by adding the following optional modules:

- **AIC-EP** for advanced control of scan patterns
- **DIA**, dynamically measure intensity over time
- **IPA**, develops icon-driven work files for automatic image analysis and processing
- **IPA-MTA**, track and analyze moving objects
- **QFA-FRET**, accurate FRET measurements and cross talk correction
- **VIS-MD**, provides rapid 3D visualization of multi-dimensional data sets
- **DNN**, Remove or Restore blur in images using fast algorithms
- **DNN-2D**, a Point Spread Function is derived and used in restoration

Getting Started

This **Quick Start Guide** contains examples of how to utilize **AIC-AD**. For further assistance, refer to the online help, manual, or visit support at <http://www.cimaging.net>, for access to the latest **How to's** and frequently asked questions. Additional support is available at e-mail: support@cimaging.net, or Tel: 412-741-7920.

Example guides:

- ◆ Add Microscope to a New Profile 2
- ◆ Add Confocal to a New Profile 3

Add Microscope to a New Profile

1. At File Menu, select **Manage Profiles** (fig.1).
 2. Click **Add** (fig.1), enter a name for the new profile.
 3. Click **Properties** to add device and file path (fig.1).
 4. Assign Default File Paths by clicking on category and browsing to desired file path (fig.2).
 5. Click on **Device Control** to add camera, microscope or confocal (fig.2).
 6. Add a camera by default clicking on **Image Capture Devices > Add > Select Camera**.
 7. Add an XY stage, filter wheel or shutter by clicking on Stage/Filter/Shutter Devices.
 8. Add a microscope by clicking on **Microscope Devices > Add** (fig.3).
 9. Select Microscope > Click on COM-Port and select Com-Port > Click on Controller > Click Auto (fig.4) to read setup details previously configured in the microscope. The components list will be updated. Click OK to accept the current settings and close the window.
- Note:** Each microscope manufacturer has its own setup software, which need to be installed and configured accurately for **AIC-AD** to configure the microscope properly.
10. Click on the **Camera** icon to Activate Capture Menu.
 11. Click on the relevant **Microscope** Tab in the Capture Window to show the components and functionality of the microscope (fig.5).
 12. The microscope is controllable from the dialogue in the window (fig.5).
 13. Select objective, cube, shutter or light source for convenient imaging; however, the microscope must be aligned for best image conditions.

Microscope components listed are accessible through the **Filter Setup** on the **Sensor** Tab to synchronize activity during automated image capture.




fig.1




fig.2




fig.3




fig.4




fig.5

Add Confocal to a New Profile

1. At File Menu, select **Manage Profiles...**
2. Click **Add (fig.1)**, enter a name for the new profile.
3. Click **Properties** to add device and assign file path.
4. Click **Properties** to add device and assign file path.
5. Assign **Default File Paths** by clicking on category, and browsing to desired file path.
6. Click on **Device Control** tab to add camera, microscope or confocal.
7. Select **Microscope Devices (fig.2)** and click **Add** Select the **VisiTech VT-Infinity Confocal Scanner (fig.3)**.
8. Click on **Auto** to automatically pull the available configurations. Check **AOTF** if an Acousto-Optic Tunable Filter (AOTF) is present (**fig.3**).
9. Click **Device List** to select a communication interface. **USB, Serial, or VT-Infinity 3**.
10. Click on **Filter Setup** and customize the **Dichroic Mirror** and **Barrier Filter** step settings (**fig.3**). Click **OK** to accept the current settings and close the window.
11. Click on the **Camera** Icon to activate **Capture Menu** window.
12. Click on the **VT-Infinity** tab in the **Capture Menu** to show the components and functionality of the confocal (**fig.4**).




fig.1




fig.2




fig.3

The VT-Infinity 3 window is summarized below:

Confocal

- Select up to four Dichroic filter positions (**fig.4.A**)
- Select up to five Barrier filter positions (**fig.4.B**)
- Home Dichroic and Barrier (**fig.4.C**)
- Select Pinhole size (**fig.4.D**)
- Conveniently Start/Stop Galvo (**fig.4.E**)
- Adjust speed of Galvo Scan Rate (**fig.4.F**)
- Adjust sweep of Galvo Scan Width (**fig.4.G**)
- Adjust Galvo sweep position with Offset (**fig.4.H**)
- Sync. Scan Width and Offset with current ROI (**fig.4.I**)

AOTF

- Adjust Laser intensity on selected lines (**fig.4.J**)
- Check On/Off up to eight Laser lines (**fig.4.K**)




fig.4

Confocal and AOTF functions are accessible through Filter Setup on the **Sensor Tab** to synchronize activity during automated image capture.