


HCIimage DIA

Getting Started Guide



Release 4.8

August 2021

This guide, as well as the software described in it, is covered under license agreement and may be used or copied only in accordance with the terms of the license agreement. The information in this manual is subject to change without notice and may not be reproduced without Hamamatsu's permission. Hamamatsu has carefully prepared this manual, however, no responsibility is assumed for possible inaccuracies or omissions. Some images are simulated.

HCIImage is a trademark of Hamamatsu Corporation. DCAM-API and ImagEM are registered trademarks of Hamamatsu Photonics K.K. (EU, Japan, U.K., U.S.A.). ORCA is a registered trademark of Hamamatsu Photonics K.K. (China, France, Germany, Japan, U.K., U.S.A.) All product and brand names are trademarks or registered trademarks of their respective companies.

Hamamatsu Corporation

360 Foothill Road, Box 6910

Bridgewater, NJ 08807-0910

USA

+1 908.231.0960

sales@hamamatsu.com

<https://camera.hamamatsu.com/us/en/index.html>

Software Support

hcsupport@hamamatsu.com

www.hcimage.com

Table of Contents

Installation

Install HCImage	3
Install DCAM-API Drivers	3

Add Devices to a Profile

Add a Camera	4
Add an Olympus IX-83 Automated Microscope	4
Add a Filter Wheel and a Shutter	6
Add a Parallel Port as an IO/LED Device	7

Calibration

Calibrate an Image from Pixels to Microns	8
Link Calibration to Objective	8
Calibrate a Stage	9

Filter Setup

Filter Wheel and Shutter Setup	10
Lambda DG-4 Filter Setup as an I/O Device	11

Capture

Capture a Color Image	13
Define a Custom SubArray for Maximum Speed	14
Control an LED using Output Trigger from the Camera	15
How to Setup a Background Subtraction	16

Sequence

Setting up a Time Lapse	17
High Speed Streaming	20

DIA Overview

Understanding the Workspace	23
-----------------------------------	----

Viewing the Data

Object Summary Statistics	27
Object Summary Graphs	28
Field Summary Statistics	29
Field Data	31

Export the Data

Export Intensity Data	32
Copy to Spreadsheet or Excel	32
Batch Export DCIMG to MPTIFF	33
Batch Export W-VIEW Images	33


DIA Analysis Examples

DIA Example	34
DIA Post Acquisition Example	35
DIA Scheduler Multi-channel Example	36
DIA Scheduler Monochrome Example	38

Troubleshooting	39
------------------------------	-----------

INSTALLATION


Install HCImage

1. Insert the HCImage installation DVD into the DVD-ROM drive. If autoplay is enabled, the HCImage setup will run automatically. If autoplay fails to start, locate your DVD-ROM drive and double-click on **setup.exe**.
2. Click **Yes**, if prompted by the User Account Controls.
3. To begin the installation wizard, click **Next**.
4. Review the Software License information and click **Yes**.
5. Review the README section for up-to-date information on software compatibility and support. When you are ready, click **Yes**.
6. On the Personalize screen, enter your registration information and click **Next**.
7. Choose the Destination Folder and click **Next**. It is recommended to install the software in the default path.
8. If you are ready to proceed with the installation, click **Install**.
9. Follow the instructions on each installation page.
10. Securely connect the dongle () to a USB port after the software installation has finished.
11. Install the appropriate DCAM-API drivers, see the instructions below, then turn your camera on prior to launching HCImage.
12. Click the **HCImage** icon on your Desktop to launch HCImage.
13. Register the software to receive technical support, please go to www.hcimage.com and click **Register**.

Install DCAM-API Drivers

Before installing the camera driver, make sure that the camera is turned off.

1. Open Windows Explorer, go to HCImage installation DVD, expand the **Drivers** folder, open the **Cameras** folder and open the **DCAM** folder. If you downloaded HCImage, please go to <http://www.dcam-api.com/> and download the DCAM-API drivers for Windows.
2. Double-click **Setup.exe** to launch the DCAM-API Setup dialog.
3. Click **Yes**, if prompted by the User Account Controls.
4. Select the appropriate driver for your Hamamatsu camera from the DCAM-API Setup dialog. If you are unsure of which driver to install, please consult the DCAM-API Compatibility Note or contact your local Hamamatsu representative. To view DCAM-API Compatibility Note, select **Others** and then click **Compatibility Note**.
5. Click **Next** to begin the installation.
6. Follow the instructions on each installation page and click **Finish** when the installation is complete.




ADD DEVICES TO A PROFILE

Add and setup peripheral hardware devices to the profile. Supported hardware includes: cameras, stages, filters, shutters, I/O and microscope devices. For a list of supported devices, please visit our website at <http://hcimage.com/support/hardware.htm>.

Add a Camera

Launch HCImage, go to File, select Current Profile and then follow the steps below to add a camera to the profile.



Add an Olympus IX-83 Automated Microscope

Olympus 64-bit Drivers from the HCImage DVD

Be advised that this application installs the Olympus Camera and Microscope drivers, as well as copying all of the necessary dlls into the HCImage directory.

1. Open the DVD contents in Windows Explorer and navigate to Drivers\Microscopes\Olympus\Olympus 3 Series\x64.
2. Double-click on **Olympus_x64 Install.exe** and follow the installation instructions.
3. Click **Yes**, if prompted by the User Account Controls

Note: If using a Hamamatsu 1394 camera, this driver installation may supersede the Hamamatsu driver causing communication problems. To recover from this issue, please see "**Unable to communicate with Hamamatsu 1394 camera**" on page 39.

Configure with the Touch Panel Controller

The microscope drivers have been installed, time to configure it using the touch panel controller (TPC) and then add it as a device in HCImage. The first step is to turn on the IX3-CBH (microscope control box) and then the touch panel controller.


Note: The "Power On" sequence for turning the equipment on before use should be: Light Source > PC > Camera > IX3-CBH > Touch Panel Controller > Launch HCImage.

An initial system setup is required when using the microscope for the first time or after replacing one of the components. The microscope is setup and configured using the TPC.

1. Go to **System Setting** in the **Menu** screen.
2. Select **Unit**, enter the components connected to the IX83 for each module and tap **OK** to save the settings.
3. Select **Optical**, enter and configure the objectives, mirror units and condenser.
4. Select **Customized**, enter the focus limits and parfocality correction.
5. When the setup is complete, tap **X** to exit to the **Menu** screen.

Add Microscope to a Profile

Once the microscope has been setup from the touch panel controller, the next step is to add the microscope to a profile and configure it in HCImage. Launch HCImage, go to File and select Current Profile. In the Device Control tab, select Microscope Devices and click Add.



The figure shows four sequential screenshots of the 'Add MICROSCOPE Device' dialog box, illustrating the steps to add a microscope device:

- 1 Device Control**: Select the microscope manufacturer and the model. The 'Controller' tab is active, showing a list of manufacturers (LEICA, NIKON, OLYMPUS, ZEISS, Hamamatsu AEQUORIA, Yokogawa CSU) and a dropdown menu set to 'IX3'. The 'Select Components' section on the right has 'Auto' selected, with checkboxes for XY Stage, Z Focus, Fluorescence Module, Lamp Control, Condenser, Dia ND Filter, Epi ND Filter, Filter Wheel No. 1, Filter Wheel No. 2, and DIA Shutter.
- 2 Device List**: Select the device that is listed. The 'Device-List' tab is active, showing a list of installed devices with 'IX3 0' selected.
- 3 Device Control**: Select Auto and HCImage will update the components list. The 'Controller' tab is active, showing the same list of manufacturers and the 'IX3' dropdown. The 'Select Components' section on the right has 'Auto' selected, with checkboxes for XY Stage, Z Focus, Fluorescence Module, Lamp Control, Condenser, Dia ND Filter, Epi ND Filter, Filter Wheel No. 1, Filter Wheel No. 2, and DIA Shutter.
- 4 Filter Setup**: Verify the settings for the individual components, modify if necessary and click OK. The 'Filter Setup' tab is active, showing a list of filter names for each installed filter position for each installed wheel. The 'FL-Block' tab is active, showing a list of filter names (Filter cube1, Filter cube2, Filter cube3) and a 'Filter Pos Limit' of 8. The 'OK' button is highlighted.

Add a Filter Wheel and a Shutter


Launch HCIImage, go to File and select Current Profile. In the Device Control tab, select Filter Devices, click Add and follow the instructions below.

The image displays four sequential screenshots of the 'Add FILTER Device' dialog box, illustrating the steps to add a filter wheel and shutter. Red dotted lines and numbered callouts (1-4) highlight the specific actions required in each step.


- 1 Device Control**
Select the controller, the type of filter wheel and enable Add Shutter.
In the first screenshot, the 'Controller' tab is active. A list of controllers is shown, with 'PRIOR Filter' selected. To the right, the 'Filter Wheel' section has 'Single' selected, and the 'Add Shutter' checkbox is checked.
- 2 Enable COM Port**
Select the COM Port for the device.
In the second screenshot, the 'COM-Port' tab is active. A list of COM ports is shown, with 'COM1' selected. The 'Communication settings' section shows 'Baud Rate: 9600', 'Data Bits: 8', 'Parity: None', and 'Stop Bits: 1'.
- 3 Filter Setup**
Select the filter position, enter the name and click Modify.
Repeat for each filter position.
In the third screenshot, the 'Filter Setup' tab is active. A list of filter positions is shown, with '380 nm' selected. The 'Customize Positions' section shows 'Wheel 1 Speed: 100' and 'Filter Pos Limit: 10'. The 'Modify' button is highlighted.
- 4 Shutter Setup**
Set the status to Active and click OK.
In the fourth screenshot, the 'Shutter Setup' tab is active. A list of shutters is shown, with 'Shutter 1' selected. The 'Shutter Settings' section shows 'Shutter 1: Active', 'Shutter 2: None', and 'Shutter 3: None'. The 'OK' button is highlighted.

Add a Parallel Port as an IO/LED Device

In the Device Control tab, select IO/LED Devices and follow the instructions below.




If the computer doesn't have a parallel port, we recommend the SIIG CyberParallel PCIe port card (<http://www.siig.com/it-products/serial-parallel/parallel/pcie/dp-cyberparallel-pcie.html>). Install the parallel port card and driver as per the instructions provided with the card and then launch the Device Manager and follow the instructions below.





CALIBRATION

Calibrate an Image from Pixels to Microns


Open or capture an image with some known distance, for example a micrometer. Click on the Calibration Properties icon ( Calibration) on the Analysis toolbar and follow the instructions below.

- 1 Add Calibration**
Click Add and enter a title (e.g., 10x)
- 2 Select Units**
Select Microns from the Units list
- 3 Calibrate the Image**
Click Calibrate. Draw a line to span the distance to measure
- 4 Enter the Distance**
Enter the known distance and click OK



Link Calibration to Objective

To link the calibration to an objective, go to the Microscope Setup panel in the Devices pane and follow the steps below.



- 1 Enable Link to Calibration**
Click and select Link to Calibration
- 2 Select Calibration**
Right-click on the 40x objective and select the 40x calibration

Calibrate a Stage

Before calibrating the stage, make sure to load the correct scale factor for the selected objective and then follow the instructions below.

1 Define Step Size
Go to the Stage Setup tab, located in the XYZ Setup panel and click Step Size

2 Center Top Left
Center an object to the green crosshair and click OK

3 Center Bottom Right
Center the same object to the green crosshair and click OK


4 Object Centered
If the object moved back to the center of the top left crosshair, the calibration was successful, click OK. If the object did not move to the same location, check the camera orientation to the microscope and repeat Steps 1-3.

FILTER SETUP

Once the filter device has been added to the profile it will need to be configured in the Filter Setup. The examples below outline the basic steps for configuring two commonly used filter devices, a filter wheel with a shutter and a Lambda DG-4.

Filter Wheel and Shutter Setup

After the filter wheel and shutter have been added to the profile, go to Filter Setup in the Device pane and follow the instructions below.



The screenshot shows the 'Filter Setup' dialog box. It has several sections: 'Enable Automated Filter/Shutter Control' with checkboxes for 'IO/LED Devices' and 'XYZ Stage Offset'; 'Advanced Settings' with checkboxes for 'Return to Idle On Exit', 'Return to Idle After Capture', 'Return to Idle During Delay', 'Dazzle Protection', and 'Exposure Protection'; 'Default Idle Positions' with a list containing 'DAPI' and buttons for 'Add', 'Copy', and 'Remove'; a 'Time Delay' section with a 'No' radio button; a 'Delay Position' dropdown set to 'Pre-Exposure'; and a table at the bottom for 'Filter-Shutter' and 'IO/LED Device' settings. A small 'Enter a Text String' dialog is open over the 'DAPI' entry in the 'Default Idle Positions' list.

- 1 Enable Return to Idle Conditions**
Select Return to Idle on Exit, After Capture and During Delay
- 2 Add Filter Group**
Click Add, enter name and click OK
- 3 Enable Filter Settings**
Right-click on the filter group that was just created and select a filter tint
- 4 Select Filter Position**
Select the filter from the list
- 5 Define Shutter Setting**
Select Open
- 6 Add Remaining Filters**
Repeat the steps to add the remaining filters
- 7 Define Default Idle Settings**
Select Default Idle Positions, under Filters select Don't Care and for Shutters select Closed

Lambda DG-4 Filter Setup as an I/O Device

TTL can be used to control many types of devices. This example explains how to configure a Lambda DG-4 as an I/O Device controlled through the parallel port. In the Device pane go to Filter Setup and follow the instruction below.

1 Enable IO/LED Devices
Select IO/LED Devices

2 Enable Return to Idle Conditions
Select Return to Idle on Exit, After Capture and During Delay

3 Add Filter Group
Click Add, enter name and click OK

4 Enable Filter Settings
Right-click the filter group that was just created and select filter tint

5 Enable IO Pin Settings
Click ellipsis button

6 Define Pin Settings
Select Don't Care (un-check), enable Pin 2 and click OK


7 Add Remaining Filters
Repeat the steps to add the remaining filters and attenuations using the pin settings in the table below

8 Define Default Idle Pin Settings
Select Default Idle Positions, then for the pin settings select Don't Care (un-check) and click OK

Filter Position	Attenuation		
	100%	50%	33%
1	Pin 2	Pins 2 & 4	Pins 2 & 5
2	Pin 3	Pins 3 & 4	Pins 3 & 5
3	Pins 2 & 3	Pins 2, 3 & 4	Pins 2, 3 & 5
4	Pin 4	Pin 5	Pins 5 & 4


CAPTURE

The Capture Pane provides a flexible and comprehensive method to access camera features and functionality. The Capture Pane is organized by functionality into panels that can be expanded when in use or collapsed when space is needed. The capture controls at the top of the pane (shown below) are always visible and used for controlling how images are acquired and displayed.



Capture Presets


Capture presets save basic settings such as the capture mode, channels, filters, exposure times, as well as output trigger settings and advanced camera properties. For a list of the camera settings that are saved, select a capture preset from the Capture Presets dialog and click View. HCImage will load the capture settings from the previous session when launched.



Note: Capture presets are not automatically saved before changing presets or exiting the software. To make changes to a saved capture preset, select the capture preset from the list, adjust the capture settings and click Save.

Capture a Color Image

Capturing a color image requires filter setup, for instructions on configuring filters, please see "Filter Setup" on page 10.




1 Select Capture Mode
Select RGB Color: 3-Band

2 Select Filters
Select Red for channel 1, Green for channel 2 and Blue for channel 3


3 Adjust Exposure
Click Live and adjust the exposure manually or use Auto Expose

4 Capture a Color Image
Click Capture1

Hint: In order to achieve the best possible speed when acquiring color images, set the same exposure for each channel. Once each of the exposures have been entered, click the Exposure Lock icon () to lock the exposure settings. Now any exposure adjustments will be made to all of the channels.

How to use AutoSave

Enabling AutoSave will automatically save the current image every time Capture1 is selected. The captured image is saved as a TIFF based on the file name and destination directory defined in the AutoSave Properties dialog. Enable AutoSave and then click on the ellipses to open the AutoSave Properties dialog.



1 TIFF or MPTIFF
Enable to save as MPTIFF for multiple image capture versus individual TIFF images

2 Set Location
Click the ellipsis icon and navigate to the destination directory

3 Set Default File Name
Enter file name

4 Save Settings
Click OK

Define a Custom SubArray for Maximum Speed

Click Live, focus on the sample and move the area of interest into the center of the image. Follow the steps below to define a custom subarray.

1

Define SubArray
Click Define and draw the area on the image

2

Center on Sensor
Right-click and select Center on Image

3

Apply SubArray
Click Apply

Binning and SubArray

Binning: 1

☒ Adjust Exposure Time

Sub-Array

Preset Sizes 2048 x 2048

Reset

Define

Apply


X0 0 Width 2048

Y0 0 Height 2048

Note: Centering the subarray for maximum speed is only required for the ORCA-Flash 4.0 series cameras.


Control an LED using Output Trigger from the Camera

Some cameras provide a range of output trigger signals to synchronize with an external instrument where the camera becomes the master and the external instrument becomes the slave.




How to Setup a Background Subtraction

Typically used in fluorescence microscopy, a background subtraction can be used when the image presents a dark non-uniform background. To perform a background subtraction click Live, bring the sample into focus and then move the stage off of the sample so that only the background is visible. Next, follow the steps below, when finished move the stage to bring the sample into view and the background subtraction is applied.




Hint: HcImage remembers the capture settings from the previous session, if background subtraction was left enabled, Process ON will be displayed in the Camera Control panel. The display image may appear distorted or black.




SEQUENCE

The Sequence pane provides a variety of options for defining a time lapse or high speed streaming. The sequence controls at the top of the pane (shown below) are always visible and used for selecting the scan type and reporting in real time, information about an ongoing sequence. This sections covers the basic steps for setting up a typical time lapse and high speed streaming.




Setting up a Time Lapse

The Scan Settings panel provides a variety of options for defining a time lapse to fit the needs of your application. This section provides three examples of typical time lapse settings, using each of the storage options.



How to Use AutoSave


In the AutoSave Properties dialog, the user can determine how and where to store the acquired data. Image data can be saved as a CXD, TIFF or MPTIFF. The example below provides a description of the Auto Save Properties dialog.



Note: MPTIFF files have a 65,000 image limit or 4 GB size limit. For image sequences exceeding these limits, multiple MPTIFF files will be saved and numbered sequentially.


Setup a Time Lapse - Save to Disk

The time lapse in this example will acquire an image every 30 seconds for 3 hours and the data will be saved as a cxd. Once you are satisfied with capture settings and the sample is in focus, go to the Sequence pane and follow the steps below.



Setup a Time Lapse - Save to the Temporary Buffer

Acquired data is stored in memory with the option to review the image sequence before saving or deleting it. When Temporary Buffer is selected, End Frame is automatically enabled and display the maximum number of frames that can be streamed to memory. Once you are satisfied with capture setting and the sample is in focus, go to the Sequence pane and follow the steps below.



1 Scan Type
Select Time Lapse

2 Auto Save
Click the ellipses icon, select CXD and enter the file location and naming convention

3 Field Delay
Select 0 Delay

4 End Frame
Enter 500

5 Temporary Buffer
Select to Temporary Buffer

6 Start Acquisition
Click Start


7 Acquisition Complete
Review acquired data using the playback controls in the Image Display

8 Save or Delete
Save - click OK
Delete - click Cancel

Note: Streaming to the Temporary Buffer is very useful because it provides the option to review the image sequence when trying to capture specific event and for demonstrating camera speeds.

Setup a Time Lapse - Save to Memory

The time lapse in this example will store images in memory until the acquisition is stopped or runs out of memory at which point the acquired images are saved to disk for the remainder of the time lapse. Once you are satisfied with capture settings and the sample is in focus, go to the Sequence pane and follow the steps below.



1 Scan Type
Select Time Lapse

2 Auto Save
Click the ellipses icon, select CXD and enter the file location and naming convention

3 Field Delay
Select 0 Delay

4 Continuous
Select Continuous


5 Memory
Select to Memory

6 Start Acquisition
Click Start

High Speed Streaming

High Speed Streaming is used to obtain the fastest acquisition speed from the camera. This scan is optimized for single channel streaming to RAM or directly to the computer's solid state drives (SSD) configured in a RAID 0.

Note: Acquisition rates will vary based on the PC configuration, for information about the computer requirements, please see the [PC Recommendations for ORCA-Flash4.0 V3 / LT+](#).



Control
Enter the number of frames to acquire and the approximate end time is displayed to the right

Stream Type
Stream directly to HDD or into memory with option to use Circular Buffer

AutoSave/AutoConvert
Define how streamed data is handled

DCIMG Location
Set a file location for streaming data to DISK

Display
Select a live display or to review acquired images

Note: High Speed Streaming does not support multi-channel acquisition, camera registration features (i.e., flip, rotation and pixel shift) or software processing operations (e.g., shade correction and rolling average).

Steps for Streaming to Disk

When streaming to disk, a temporary file (.dcimg) is created to store the data while it is being acquired, the temporary file location needs to be located on the RAID array, SSD drive, or the fastest drive available. Configure the capture settings, go to the Sequence pane and follow the steps below.

The screenshot shows a software interface for configuring streaming to disk. It includes a 'Select Scan Type' dropdown set to 'High Speed Streaming', a 'Start' button, and a 'Progress' section with a frame counter and time elapsed. The 'Scan Settings' section is expanded, showing 'Frame Count' set to 1000, 'Best Time' at 9.9003 sec, and 'DISK' selected as the stream type. The file path is 'D:\Experiment Data\DCIMG\rec*.dcimg'. Below this, 'AutoConvert' is checked, and file format options (CXD, TIFF, MPTIFF) and output modes (Live Image, Review) are visible. Five numbered callouts provide step-by-step instructions: 1. Select Scan Type (High Speed Streaming), 2. Enter Frame Count (1000), 3. Select Stream Type (DISK), 4. Auto Convert File Type (checked), and 5. Start Streaming (click Start).

- 1 Select Scan Type**
Select High Speed Streaming
- 2 Enter Frame Count**
Enter the number of images to acquire
- 3 Select Stream Type**
Select DISK
- 4 Auto Convert File Type**
Enable AutoConvert and select file type
- 5 Start Streaming**
Click Start

Note: To leave the streamed data as a DCIMG file disable AutoConvert.

Steps for Streaming to RAM

Acquired data is stored in memory with the option to review the image sequence before saving or deleting it. In the AutoSave Properties dialog, the user can determine how and where to store the acquired data. Once you are satisfied with capture settings and the sample is in focus, go to the Sequence pane and follow the steps below.

Note: The Circular Buffer stores streamed data in memory, once the frame count has been reached, the previous acquired data is replaced sequentially. The cyclic process repeats until the acquisition is stopped, leaving the most recent images stored in RAM.

The screenshot shows the 'AutoSave Properties' dialog box with the following settings and annotations:

- 1 Select Scan Type:** A red box highlights the 'High Speed Streaming' dropdown menu.
- 2 Enter Frame Count:** A red box highlights the 'Frame Count' input field, which contains the value '1000'.
- 3 Select Stream Type:** A red box highlights the 'RAM' radio button under the 'Stream Type' section.
- 4 Auto Save File Type:** A red box highlights the 'AutoSave' checkbox and the file type options (CXD, TIFF, MPTIFF).
- 5 Start Streaming:** A red box highlights the 'Start' button.


Other visible settings include: 'Best Time' at 9.9003 sec, 'Event Marker' set to 0, and 'Delay Remaining' at 00:00:00.

DIA OVERVIEW

Dynamic Intensity Analysis (DIA) is optimized for high speed processing and intensity analysis over time, including Live viewing of images and data simultaneously. Measuring and plotting of data is available on-line or off-line, and may be access by clicking DIA Analysis in the Sequence Pane. This functionality is only available in HCImage DIA and HCImage Analysis.


Understanding the Workspace

The Side Panel includes the Dynamic Intensity Analysis functionality, that is accessed through the Sequence pane by selecting DIA Analysis. Once enabled, the Intensity Analysis and Graph Setup panels are available, providing the tools to setup an experiment without having to switch panes.




Scan Settings

The Scan Settings panel is easy-to-use, simply set the speed, define the capture interval, enter the number of images to capture and where to save the data.




Intensity Analysis

The new Intensity Analysis panel is configured based on the selected Analysis mode: Simple or Advanced. The Advanced mode provides a comprehensive set of tools to help identify large numbers of objects and objects that are not easily differentiated. The Simple mode provides a variety of drawing tools that can be used to manually identify objects of interest.



To switch between the two analysis modes go to **View** on the menu bar, then highlight **Analysis Mode** and select **Advanced**.




In the Advanced Mode, the user has the ability to identify a large number of objects, as well as, identify hard to detect objects because of defects in the image. The tools are grouped by function as described below:

- **Enhance** the image for detection by correcting for defects such as noise, dust, intensity variation, etc.
- **Identify** objects using an intensity threshold and create a binary image layer over the image.
- **Modify** the binary image layer, filling holes and separating connected objects.
- **Qualify** objects based on one or more measurement ranges and conditionally remove edge objects.
- **Draw/Erase** objects using a set of drawing tools.

Measure Objects

Intensity measurements are available in the Select Intensity Measurements dialog, click the Select Measurements icon to open the dialog. Select measurements by clicking the measurement check box to the left. Filter the view of the measurements by selecting one of the categories in the right. When correcting for background fluorescence, the corrected and uncorrected data for each of the selected measurements will be collected.



Note: When using the ORCA-Flash4.0 LT, ORCA-Flash4.0 V2/V3 or the ImagEM X2, select Measure in Electrons, to report the intensity measurement values in electrons. When Electrons is selected, measured data will ONLY be reported in electrons. Measurement names will be preceded by an "e" denoting the measurement is in electrons.

Custom measurements are available to deal with complex situations, use the built-in equation editor to apply standard measurements and mathematical functions to customize special measurements to suit specific applications. To create a custom measurement, click Select Measure and then click Custom to open the equation editor.

Graph Setup

The Graph Setup panel lets users decide the measurements that will be displayed during the experiment and how the graphs will be displayed. The user can choose to display the data from a single object, the average of all of the objects, or all of the objects. The data for all of the measurements in the View Measurements list will be collected and saved, regardless of whether they are graphed during the experiment. Also, keep in mind that depending on the number of objects and measurements selected, the graph will become very crowded and it may become hard to differentiate the objects.

Image
Select which image to display in the graph when multiple monochrome images are acquired

Object
Select which objects will be displayed on the graph

Group List
A list of objects organized by size and class

Graphs
Select which graph settings to display

Class
Select which class to display or choose all classes

Group Size
Enter the number of objects to be included in a group

Display History
Select whether to display the entire graph or a specific number of fields or time segment. Choose Fit, to resize the graph to fit the window rather than expand in blocks.

Image #: 1

Class #: 2

Object: Group List

Group: 2.1 - 2.5

Size: 5

Count: 15 <-Count

Selected: 5 <-Selected

View Measurements:

- ☒ MEAN_RED
- ☒ MEAN_GREEN
- ☒ RATIO_OF_MEANS_RG

☒ Raw Data

☐ Corrected Data

Display History:

- ☒ Full
- ☐ Partial
- ☐ Fit

100 Fields 30 sec

- ☒ X Axis = Field Number
- ☐ X Axis = Field Time

Edit Group
Select which objects to display from a list of objects in a group

Group Average
Display the mean value of the selected group of objects

View Measurements
Select which measurements to plot on the graph and which data to display:
Raw - uncorrected object data
Corrected - raw object data with the background fluorescence signal removed


Show Legend and Events
Select whether to display the Legend and show Events during the experiment. Keep in mind that as the number of objects and measurements increase, the legend will become very large.

X Axis
Display the X axis as the Field Number or the Field Time

VIEWING THE DATA

Object Summary Statistics

Object Summary Statistics are collected for each of the Object Measurements made for each Measurement Class. As Object Measurements in a Workfile may be selected and deselected during data collection the Count value may vary between Field Measurements. Each statistic is computed according to the actual count of objects processed for each measurement selected.



STATISTIC	MEAN_RED.1	MEAN_RED_corrected.1	MEAN_GREEN.1	MEAN_GREEN_corrected.1	RATIO_OF_MEANS
Minimum	6423.400000	1435.312009	13740.924712	56.832770	0.241989
Maximum	13941.772682	8790.932979	36978.963636	23988.617957	0.780873
Mean	9528.211858	4523.514024	20550.985124	7110.994569	0.493426
Smp Std Dev	1750.612391	1773.226967	5123.465532	5118.109874	0.150743
Total	41847906.482077	19867273.593188	90259926.666736	31231488.148217	2167.128090
Smp Variance	3064643.744883	3144333.878158	26249899.059030	26195048.677956	0.022723
Pop Std Dev	1750.413084	1773.025086	5122.882226	5117.527178	0.150726
Pop Variance	3063945.966253	3143617.955144	26243922.306056	26189084.413685	0.022718
Std Error	26.415500	26.756738	77.309462	77.228649	0.002275
Mean Variance	697.778630	715.923014	5976.752973	5964.264271	0.000005
Sqr Total	412192569475.373110	103676660775.431240	1970193717033.3...	337109401356.980470	1169.096482
Recip Total	0.476154	1.149833	0.225817	1.203564	9858.296827
Count	4392	4392	4392	4392	4392


Statistical Measurements

The Statistics computed for Object and Field Measurements are derived as follows:

Statistic	Note	Formula
Count of Items	total number of items considered in the data set	Count = N
Minimum	minimum occurring value in the data set	Min (x)
Maximum	maximum occurring value in the data set	Max (x)
Total Value	sum of all values in the data set	Total value = ($\sum x$)
Mean	total value divided by count of items	$\mu = (\sum x)/N$
Sample Variance	used to characterize incomplete samples	$s^2 = (\sum x^2 - (\sum x)^2/N)/(N-1)$
Sample Standard Deviation	used to characterize incomplete samples	$s = \sqrt{s^2}$
Population Variance	used to characterize complete samples	$\sigma^2 = (\sum x^2 - (\sum x)^2/N)/N$
Population Standard Deviation	used to characterize complete samples	$\sigma = \sqrt{\sigma^2}$
Standard Error of the Mean	experimental uncertainty of an averaged measurement	$SE_{\mu} = s / \sqrt{N}$
Total of Values Squared	sum of squares	$\sum x^2$
Total of Reciprocal Values	sum of reciprocals	$\sum 1/x$

Object Summary Graphs

Object Measurements can be plotted for each object in the data document. The Object Summary Graphs show Object Measurement data of all fields. The interactive graph lets the user customize the display. Use the graph toolbar to select the measurements and which class and objects to display. Use the right-click menu to customize the look and feel of the graph by changing the title and legend fonts or adjusting background and border colors. The right-click menu also allows users to display the X axis as number of fields or field time.




Object Summary Graph Toolbar


The toolbar provides multiple options for managing how the data is displayed.

<p>Image Select which image to display in the graph when multiple monochrome images were acquired</p>	<p>Object Select which objects will be displayed on the graph</p>	<p>Edit Group Select which objects to display from a list of objects in a group</p>	<p>Group Size Select or enter the number of objects to be included in a group</p>	<p>Data in Y Axis Display the measurement data in the Y axis</p>
<p>Event Markers Display event markers</p>	<p>Class Select which class to display or choose all classes</p>	<p>Group List A list of objects organized by size and class</p>	<p>Group Average Display the mean value of the selected group of objects</p>	<p>Show Legend Display legend</p> <p>Select Measurements Select which measurements to display in the graph</p>

Object Summary Data

In addition to the Object Summary graphs, the intensity measurement data can also be displayed using a Spreadsheet View and a Table View. To change the view, go to the Image Data Views


toolbar, click on the Current View icon () and select either Spreadsheet View or Table View.



The left screenshot shows the 'Gluc 04201622.cxd' window with the 'Object Summary Statistics' node selected. The data is displayed in a spreadsheet format with columns for FIELD#, EVENT MARKER, FLD TIME (sec.), MEAN_RED.1, MEAN_GREEN.1, RATIO_OF_MEANS_RG.1, and MEAN_OF_RATIOS_RG.1. The right screenshot shows the same data in a table format, with the same columns and data values.

Field Summary Statistics

Field Summary Statistics are collected for each of the Field Measurements made for each Measurement Class. As Field Measurements in a Workfile may be selected and deselected during data collection the Count value may vary between Field Measurements. Each Statistic is computed according to the actual Count of Fields processed for each measurement selected.




The screenshot shows the 'Field Summary Statistics' node selected in the 'Data Tree' for 'Gluc 04201621.cxd'. The statistics are displayed in a table format with columns for STATISTIC, MEAN_RED.1, MEAN_RED_corrected.1, MEAN_GREEN.1, MEAN_GREEN_corrected.1, and RATIO_OF_MEANS_1.

STATISTIC	MEAN_RED.1	MEAN_RED_corrected.1	MEAN_GREEN.1	MEAN_GREEN_corrected.1	RATIO_OF_MEANS_1
IMAGE 1					
Minimum	8412.867153	3439.850790	16527.178023	2825.896642	0.326850
Maximum	11267.487605	6125.431425	25910.451251	12295.486170	0.678166
Mean	9975.108938	4970.411103	21013.210314	7573.219759	0.496595
Smp Std Dev	1074.311789	1109.583456	3598.762258	3558.240672	0.130058
Total	1825444.935636	909585.231932	3845417.487455	1385899.215850	90.876802
Smp Variance	1154145.820494	1231175.445874	12951089.788437	12661076.677601	0.016915
Pop Std Dev	1071.372490	1106.547654	3588.916105	3548.505385	0.129702
Pop Variance	1147839.012732	1224447.711197	12880318.805987	12591890.466248	0.016823
Std Error	79.415413	82.022769	266.028161	263.032719	0.009614
Mean Variance	6306.807762	6727.734677	70770.982450	69186.211353	0.000092
Sqr Total	18419066632.441399	4745086467.485528	83161664750.384...	12800035280.429529	48.207483
Recip Total	0.018570	0.038972	0.008963	0.031482	397.220220
Count	183	183	183	183	183


Opening the Field Summary Statistics node will display a node for each Measurement Class present. Selecting the Field Summary Statistics node will display all Field Measurement Classes. Under the Field Summary Statistics node is a node for each Measurement Class. Selecting the Class node will display the Field Measurement Data for the individual Class.


Field Summary Graphs

Field Measurements can be plotted for each Field in the Data Document. The Field Summary Graphs show Field Measurement parameters for each Measurement Class with the Field number as the X axis variable and the Field Measurement as the Y axis variable.



The Field Summary graphs are an easy to use interactive way of displaying and viewing data. First the Mean Red (Corrected 340nm) intensity is measured over time. Second there are two Event Markers that indicate specific points during the experiment that an event happened. In this case, we know the name of the reagents that were added and the time they were added. The Event

Markers may be toggled on/off by clicking the Show Event Markers Icon (). The intensity at any given point is displayed by hovering the cursor over a section of the graph. The corresponding image may also be viewed by clicking on a point along the graph. Zoom in on a specific area of the graph by clicking and dragging the mouse over the area of interest. Release the button and the graph zooms to the size of the box. Click the left mouse button to return to the normal view. While in a zoomed position the user can pan by dragging the mouse in the direction they wish to view. Using the features in the Playback Toolbar we can play the image sequence and visually see the changes in the intensity and how they are plotted on the corresponding graph. Use the right-click menu to customize the look and feel of the graph by changing the title and legend fonts or adjusting background and border colors. The right-click menu also allows users to display the X axis as number of fields or field time.



Field Data

Field Data contains information for each field about when and where the image was captured. The data can be viewed in a Table View or Spreadsheet View and copied to the Windows Clipboard. These details can include:

- X,Y,Z Stage Position Microns
- Image Width in pixels
- Image Height in pixels
- Image Depth in bits per pixel
- Time From the Start (Hours:Minutes:Seconds.Hundredths)
- Time From Last (Hours:Minutes:Seconds.Hundredths)
- Computer-controlled Wavelength used (in nanometers)
- Group Number
- Group Index
- Group Size

Gluc 04201621.cxd

Fld#...	FldName	Event_Marker	Time_From_Start	Time_From_Last	MEAN_RED.1	MEAN_RED_corrected.1
53			0:01:23.890999	0:00:1.625000	8640.913408	3566.641803
54		lonmycin	0:01:25.500000	0:00:1.609000	8568.640913	3439.850790
55			0:01:27.110001	0:00:1.610000	8652.040742	3481.719755
56			0:01:28.719002	0:00:1.609000	8713.363011	3556.054369
57			0:01:30.328003	0:00:1.609000	8774.581986	3568.581986
58			0:01:31.985001	0:00:1.657000	8892.187115	3680.755016
59			0:01:33.610001	0:00:1.625000	9032.991118	3714.324451
60			0:01:35.203003	0:00:1.593000	9126.335076	3934.878286
61			0:01:36.813004	0:00:1.610000	9223.989395	4036.878284
62			0:01:38.438004	0:00:1.625000	9288.091550	4069.819946
63			0:01:40.046997	0:00:1.609000	9210.661075	4092.340087
64			0:01:41.656998	0:00:1.610000	9239.303845	4136.414956
65			0:01:43.265999	0:00:1.609000	9290.299666	4168.003370
66			0:01:44.875000	0:00:1.609000	9351.068198	4221.907705


Blob Count: 24

Zoom: 100% SF: 1.000Px Default

EXPORT THE DATA


Export Intensity Data

Use this method to export all of the collected intensity data from the data document to a spreadsheet. This includes the object and field data as well as the object and field summary statistics. With the data document open follow the steps below to export the intensity data to a spreadsheet.




Copy to Spreadsheet or Excel

To copy only the data from a specific graph, select the graph and follow the instructions below.



Batch Export DCIMG to MPTIFF

In the File menu select Batch Export and follow the instructions below. The exported files are not automatically opened in the software.




- 1 Enter Source Location**
Type: Select DCIMG Files
Browse: Go to the file directory
- 2 Enter Destination Location**
Type: Select Multi-Page TIFF Files
Browse: Go to output directory
- 3 Define Output File Name**
Define the file naming convention
- 4 Enable Create Series Folder**
Select Create folder for TIFF series
- 5 Export to MPTIFF**
Click OK

Note: MPTIFF files have a 65,000 image limit and 4 GB size limit. For image sequences having more than 65,000 images or larger than 4 GB, multiple MPTIFF files will be saved and numbered sequentially.

Batch Export W-VIEW Images

In the File menu select Batch Export and follow the instructions below.



- 1 Enter Source Type and Location**
Type: Data Select Data Documents
Browse: Go to the file directory
- 2 Enter Destination Type and Location**
Type: Select TIFF Files
Browse: Go to the output directory
- 3 Enable Create Series Folder**
Select Create folder for TIFF series
- 4 Define Channel Options**
Enable Split Image and select Single Color Image, A-Red B-Green and Top/Bottom
- 5 Export to TIFF**
Click OK

DIA ANALYSIS EXAMPLES

DIA Analysis can be run live, in real time, as well as on previously acquired data sets. The following examples provide instructions for various situations.

DIA Example

The instructions below outline the steps for setting up a basic DIA Analysis experiment with a single class of objects and background correction. Configure the capture settings as needed for your sample. Two channel, Red and Green settings were used for this example.

The screenshot shows the DIA Analysis software interface with the following components and steps:

- 1 Enable DIA Analysis**
Select DIA Analysis (checked).
- 2 Define Scan Settings**
Expand the Scan Settings panel, select Field Delay 1 and enter 5 s.
- 3 Identify Objects of Interest**
Expand the Intensity Analysis panel, click on the ellipse icon and identify the cells of interest in the Image Display.
- 4 Identify the Background**
Select Background, click on the ellipse icon and identify a background area in the image.
- 5 Define the Measurements**
Click on the Select Measurements icon, select Mean Red, Mean Green, Ratio of the Means RG and click OK.
- 6 Configure Graph 1 Settings**
Expand the Graph Setup panel and disable Ratio of Means RG.
- 7 Configure Graph 2 Settings**
Select Graph 2, disable Mean Red, Mean Green and enable Corrected Data.
- 8 Start Acquisition**
Click Start

The interface includes the following sections:

- Select Scan Type:** Time Lapse (dropdown), Use Scheduler (unchecked), DIA Analysis (checked).
- Progress:** 0 fps, Time Elapsed: 00:00:00.00, Delay Remaining: 00:00:00.
- Event Marker:** 0 (selected), 1, 2, 3, 4, 5, 6, 7, 8, 9.
- Scan Settings:** Expandable panel.
- Intensity Analysis:** Object Classes (Class 1, Background 1), Add, Background, Remove, Ratio, Select Measure, Show All Blobs (unchecked).
- Draw Objects:** Draw (selected), Erase, and various drawing tools.
- Graph Setup:** Graph 1 (selected), Graph 2 (selected), Visible (checked), Image #: 1, Object: Average, Class #: 1, Group (selected), Size: 10, View Measurements: MEAN_RED (unchecked), MEAN_GREEN (unchecked), RATIO_OF_MEANS_RG (checked), Raw Data (checked), Corrected Data (checked), Display History: Full (selected), Partial (30 seconds), Show Legend (checked), Show Events (checked), X Axis = Field Time (selected).

DIA Post Acquisition Example

DIA Analysis can be run on previously acquired image sequences. The example below includes a single class of objects with background subtraction. Open the data document to analyze, go to the Sequence pane and follow the instructions below.


Note: In addition to data documents (cxd), DIA Analysis can run directly from multi-page tiff and dcimg files.

The screenshot shows the HCImage DIA Analysis interface. The interface is divided into several panels: 'Select Scan Type', 'Progress', 'Intensity Analysis', and 'Graph Setup'. Red dotted lines with numbered circles (1-9) point to specific controls and dialog boxes, indicating the sequence of steps for running the analysis.

- 1 Enable DIA Analysis and Select Data Set**
Select DIA Analysis and then select the data set from the Select Scan Type list.
- 2 Identify Objects of Interest**
Expand the Intensity Analysis panel, click on the ellipse icon and identify the cells of interest in the data set.
- 3 Identify the Background**
Select Background, click on the ellipse icon and identify a background area in the image.
- 4 Define the Measurements**
Click on the Select Measurements icon, select Mean Red, Mean Green, Ratio of the Means RG and click OK.
- 5 Configure Graph 1 Settings**
Expand the Graph Setup panel and disable Ratio of Means RG.
- 6 Configure Graph 2 Settings**
Select Graph 2, disable Mean Red, Mean Green and enable Corrected Data.
- 7 Start Acquisition**
Click Start
- 8**
A dialog box titled 'HCImage' appears, asking to 'Select the output format for your analysis.' The 'DataDoc' option is selected.
- 9**
A dialog box titled 'HCImage' appears, asking 'Do you want to replace data in "D:\Fura gluc.cxd" or create a new file?'. The 'Replace' option is selected.

DIA Scheduler Multi-channel Example

The diagram below provides a partial timeline of the experiment. Each component is broken down and plotted as interval time (red line), with the individual intervals labeled by action. Components are grouped into scans (purple line) and labeled by type. This example provides instructions for creating the schedule and defining the scans, identifying objects of interest, selecting measurement and graph setup have already been covered and are not included in the DIA Scheduler example.




The schedule for the experiment consists of four scans, the scans are defined below:

1. **Scan 0** (Acquire) 340 nm and 380 nm with 50 ms exposures and a field delay of 900 ms for 15 minutes
2. **Scan 1** (Pulse) excite at 500 nm for 350 ms single pulse
3. **Scan 2** (Wait) delay of 650 ms
4. **Scan 3** (Acquire) 340 nm and 380 nm with 50 ms exposures and a field delay of 900 ms for 10 minutes

Create the Schedule

Before setting up the scans we need to enter the capture settings and configure the TTL settings for triggering pulses. For this example, the light source is controlled using TTL, where pin 2 (340 nm), pin 3 (380 nm), and pin 4 (500 nm).

Capture Settings



1

Define Capture Settings

Enable RGB Color: 2-Band, select 340 nm and 380 nm filters and enter a 50 ms exposure for each channel.

I/O Setup

Click the **I/O Setup** button in the Devices pane, the TTL trigger outputs are identified as Scan Pulse On, 0-3. If the I/O Setup button is not visible you will need to add an I/O device in the Profile.

Input and Output Job Settings

Link each I/O pin to a specific Capture Event or an Event Marker key. In Event Markers tab, you can customize names by highlighting the individual and press F2 key.

Enable I/O During Capture
☒ Capture Trigger/Signal
☐ Event Marker

Capture Events			
Capture Inputs	Input Pins	State	Frequency
Sequence Start	None	None	None
Sequence Stop	None	None	None
0: Scan Start	None	None	None
1: Scan Start	None	None	None
2: Scan Start	None	None	None
3: Scan Start	None	None	None
4: Scan Start	None	None	None

Capture Outputs	Output Pins	State
Sequence On	None	None
Field On	None	None
Exposure On	None	None
0: Scan Pulse On	LPT1 P4-	High
1: Scan Pulse On	None	None
2: Scan Pulse On	None	None
3: Scan Pulse On	None	None

OK Cancel

- 2 Define TTL Output**
 For 0: Scan Pulse On, click None under Output Pins and select LPT1 P4.

Create Scans

Go to the Sequence pane, select Use Scheduler and DIA Analysis and follow the instructions below.

Select Scan Type: Time Lapse ☒ Use Scheduler ☒ DIA Analysis

Frame Bundle: 2 Start Stop

Progress: 0 fps Time Elapsed: 00:00:00.00

Event Marker: 0 Delay Remaining: 00:00:00

Schedule Loop: 1

Name	Type	Duration	Progress	Status
Scan 0	Time	900.00 s	% 0	Not Started
Scan 1	TTL Out	0.35 s	% 0	Not Started
Scan 2	Delay	0.65 s	% 0	Not Started
Scan 3	Time	900.00 s	% 0	Not Started

1 Define Scan 0 - Acquire
 Double-click Time and enter the settings below. Right-click and select Add Subscan.

Edit SubScan

Setup
 Field Delay: 900.0 ms
☐ End Frame: 1
☒ End Time: 15.0 min

☒ to Disk
☐ to Memory (214) RAM...

OK Cancel

2 Define Scan 1 - Pulse
 Select TTL Out from the list, double-click TTL Out and enter the settings below. Add the next Subscan.

Edit SubScan

TTL Out: 0 ☐ Asynchronous
☐ Wait: 0.0 sec before 1st Pulse

Setup
 Pulse Interval: 0.0 sec
☒ End Count: 1
 Pulse Duration: 350.0 ms
☐ End Time: 0.0 sec

OK Cancel

- 4 Define Scan 3 - Acquire**
 Double-click Time and enter the settings below.

Edit SubScan

Setup
 Field Delay: 900.0 ms
☐ End Frame: 1
☒ End Time: 10.0 min

☒ to Disk
☐ to Memory (261) RAM...

OK Cancel

- 3 Define Scan 2 - Wait**
 Select Delay from the list, double-click Delay and enter the settings below. Add the next Subscan.

Edit SubScan

Setup
 Delay: 650.0 ms
☐ End Frame: 1
☐ End Time: 0.0 sec

OK Cancel

DIA Scheduler Monochrome Example

This example provides instructions for defining the scans, identifying objects of interest, selecting measurement and graph setup have already been covered and are not included in the DIA Scheduler example.

Create Scans

Go to the Sequence pane, select Use Scheduler and DIA Analysis and follow the instructions below.

The screenshot shows the DIA Scheduler interface with the following components:

- Select Scan Type:** Time Lapse (dropdown), ☒ Use Scheduler, ☒ DIA Analysis.
- Frame Bundle:** 2 (dropdown), Start, Stop buttons.
- Progress:** 0 (input), fps (input), Time Elapsed: 00:00:00.00.
- Event Marker:** 0 (dropdown), Delay Remaining: 00:00:00 (input).
- Schedule:** Loop: 5 (dropdown), table with 5 columns: Name, Type, Duration, Progress, Status.
- Scan Settings:** ☒ AutoSave, ☒ CXD, ☐ TIFF, ☐ MPTIFF, ☒ Live Image, ☐ Review.

Five numbered callouts provide instructions:

- 1 Define Loops**
Set the numbers of Loops to five.
- 2 Define Scan 0 - Acquire**
Double-click Time and enter the settings below. Right-click and select Add Subscan.
- 3 Define Scan 1 - Wait**
Select Delay from the list, double-click Delay and enter the settings below. Add the next Subscan.
- 4 Define Scan 2 - Acquire**
Double-click Time and enter the settings below.
- 5 Define Scan Settings**
Select Average each Loop.


The **Edit SubScan** dialog boxes show the following settings:


- Scan 0 - Acquire:** Field Delay: 5 sec, ☐ End Frame: 1, ☒ End Time: 5 min, ☒ to Disk, ☐ to Memory (685), RAM...
- Scan 1 - Wait:** Delay: 1 min, ☒ End Frame: 1, ☐ End Time: 0.0 sec.
- Scan 2 - Acquire:** Field Delay: 0.0 sec, ☐ End Frame: 1, ☒ End Time: 1 min, ☐ to Disk, ☒ to Memory (685), RAM...

TROUBLESHOOTING

System Information

What version of Windows is installed?


Press the **Windows Logo**  + **Pause/Break** keys to view the System Properties window. The Windows System Properties displays a basic overview of the computer including Windows edition and System type (i.e., 32-bit or 64-bit).




Unable to communicate with Hamamatsu 1394 camera

Was the Hamamatsu 1394 driver re-installed after installing the Olympus drivers?


1. Right-click on **My Computer**, select **Manage** and select **Device Manager** in the System Tools list.
2. Right-click on the appropriate 1394 controller and select **Update Driver Software**.



3. Choose **Browse my computer for driver software**
4. Choose **Let me pick from a list of device drivers on my computer**



- Choose the DCAM compatible driver - **LSI 1394 OHCI Compliant Host Controller** and click **Next** to install the driver.



Light remains on during delay or after capture

Is the light source shuttered or is it an IO/LED device?

Review the Filter Setup setting below. The example on the left is for a shuttered light source and the example on the right is for an IO/LED device.

