

# HCIimage

## CoolLED pE-4000 Universal Illumination System



© 2017 Hamamatsu Corporation.  
All rights reserved.

This guide, as well as the software described in it, is covered under license agreement and may be used or copied only in accordance with the terms of the license agreement. The information in this manual is subject to change without notice and may not be reproduced without Hamamatsu's permission. Hamamatsu has carefully prepared this manual, however, no responsibility is assumed for possible inaccuracies or omissions. Some images are simulated.

HCImage is a registered trademark of Hamamatsu Corporation. DCAM-API and ORCA are trademarks of Hamamatsu Photonics K.K. All product and brand names are trademarks or registered trademarks of their respective companies.

***Hamamatsu Corporation***

*360 Foothill Road, Box 6910*

*Bridgewater, NJ 08807-0910*

*USA*

*+1 908.231.0960*

*sales@hamamatsu.com*

*http://sales.hamamatsu.com/*

***Software Support***

*hcsupport@hamamatsu.com*

*www.hcimage.com*

# Table of Contents

---

## **CoolLED pe-4000**


Wavelength Grouping .....	4
CoolLED Driver .....	5
Add pE-4000 to a Profile .....	5
Filter Setup Examples .....	7

## COOLED PE-4000

The CoolLED pE-4000 has 16 selectable LED wavelengths ranging from 365 nm to 770 nm, that can be matched to the filters and fluorophores of almost any microscope. Precise intensity control of selected LED wavelengths in 1% steps (0–100%). The pE-4000 does not require shutters, the LEDs provide instant on/off and there isn't a warm up or cool down period.

### Wavelength Grouping

The pE-4000 provides 16 selectable wavelengths arranged in 4 channel groups, which allows for matching with available dual, triple and quad filter sets. Each channel can be individually controlled, both in terms of intensity and wavelength switching. This allows for up to a simultaneous 4 wavelength excitation, one wavelength per channel.



pE-4000 Selectable Wavelengths


Channel A	Channel B	Channel C	Channel D
365 nm	460 nm	525 nm	635 nm
385 nm	470 nm	550 nm	660 nm
405 nm	490 nm	580 nm	740 nm
435 nm	500 nm	595 nm	770 nm

**Note:** The LEDs are grouped so that any chosen combination of four will have a wavelength selected from one of 4 separate groups in order to prevent the excitation and emission spectra from overlapping.

## CoolLED Driver

For USB operation, the CoolLED Driver is required for setting up the 2 virtual COM ports used by the pE-4000. Copy and unzip the CoolLED-pE-inf.zip to a folder on the desktop. Contact CoolLED if you do not have a copy of the driver.


1. Turn on the pE-4000, let it initialize and then connect the USB cable to the computer.
2. The Windows Driver installation will fail, go to the Device Manager, right-click on **My Computer**, select **Manage** and select **Device Manager** in the System Tools list.
3. Go to **Other Devices**, right-click on **USB Virtual Serial Port A**, select **Update Driver Software...** and choose **Browse my computer for driver software**.
4. Point to the **CoolLED-pE-inf** folder on the desktop, click **Next** and allow Windows to install the driver.
5. The pE-4000 is now listed under **Ports (COM & LPT)** as **CoolLED pE-4000 USB Virtual Serial Port A (COM#)** with an assigned COM port number.
6. Go to **Other Devices**, right-click on **USB Virtual Serial Port B** and repeat the steps to update the driver for this COM port.




## Add pE-4000 to a Profile

The next step is to add the pE-4000 as an IO/LED Device to a profile in HCImage.

1. Launch HCImage, go to **File** and select **Current Profile**. In the **Device Control** tab, select **IO/LED Devices** and click **Add**.
2. Select **CoolLED precisExcite** from the list, go to the **COM-Port** tab and select the COM # assigned by Windows.



- Click **OK** to return to the Current Profile. The pE-4000 is now listed as a the CoolLED-pE under IO/LED Devices. Click **OK** to save the device settings to the profile and close the window.



- Go to the **Devices** pane and expand the **Filter Setup** panel to define filter groups. Click on the **IO/LED Device** tab to access the filter and shutter controls for the pE-4000.

Filter-Shutter	IO/LED Device
CoolLED-pE Filter-B	Don't care
CoolLED-pE B	Don't care
CoolLED-pE Filter-C	550
CoolLED-pE C	20.0 %
CoolLED-pE Filter-D	Don't care
CoolLED-pE D	Don't care
<b>Shutters</b>	
CoolLED-pE Shutter-A	Low
CoolLED-pE Shutter-B	Low
CoolLED-pE Shutter-C	High
CoolLED-pE Shutter-D	Low

**Channel Wavelength**  
Each channel has 4 wavelengths to choose from

**Channel Intensity**  
Set the intensity 0-100% (1% steps) for the selected wavelength

**Shutter Controls**  
Set the shutter state to Low (Closed), High (Open) or Don't Care (Ignore)

## Filter Setup Examples

Follow the steps below to configure the channel settings for pE-4000.

### Part 1 - Filter Setup for Default Idle Positions


The screenshot shows the 'Filter Setup' dialog box with the following settings and callouts:

- 1 Automated Control**:  Enable Automated Filter/Shutter Control. Enable IO/LED Device control.
- 2 Advanced Settings**:  Return to Idle On Exit,  Return to Idle After Capture,  Return to Idle During Delay. Enable Return to Idle On Exit, Return to Idle After Capture and Return to Idle During Delay.
- 3 Filter Group**:  Default Idle Positions. Select Default Idle Positions.
- 4 IO/LED Device Controls**:  IO/LED Device. Select to display IO/LED Device controls.
- 5 Shutter Controls**: CoolLED-pE Shutter-A, CoolLED-pE Shutter-B, CoolLED-pE Shutter-C, CoolLED-pE Shutter-D. Set the shutter state to Low for Shutters A-D.


Time Delay:  None,  Manual,  Automatic: 0.1 Sec. Delay Position: Pre-Exposure

Filter-Shutter	IO/LED Device
CoolLED-pE Filter-C	Don't care
CoolLED-pE C	Don't care
CoolLED-pE Filter-D	Don't care
CoolLED-pE D	Don't care
<b>Shutters</b>	
CoolLED-pE Shutter-A	Low
CoolLED-pE Shutter-B	Low
CoolLED-pE Shutter-C	Low
CoolLED-pE Shutter-D	Low


## Part 2 - Filter Setup for 550 nm



- 1 Filter Group**  
Select Default Idle Positions
- 2 Copy Filter Group**  
Click Copy, enter 550 nm and click OK
- 3 Filter Group**  
Right-click on 550 nm, select Customize from the tint list and follow the instructions for **A** below to create a custom tint
- 4 Filter Controls**  
Set Filter C to 550, select C, click the ellipses and follow the instructions for **B** below to enter the intensity value
- 5 Shutter Controls**  
Set Shutter C to High (Open)



- 1 Add Custom Tint**  
Click Add, enter the tint name and click OK
- 2 Tint Wavelength**  
Enter the wavelength
- 3 Save Settings**  
Click OK



- 1 Enable Voltage Settings**  
Deselect Don't Care
- 2 Intensity**  
Enter the intensity 0-100%
- 3 Save Settings**  
Click OK



## Part 3 - Filter Setup for White Light

The screenshot shows the 'Filter Setup' window with the following settings:

- Enable Automated Filter/Shutter Control:**
- IO/LED Devices:**  **XYZ Stage Offset:**
- Advanced Settings:**
  - Return to Idle On Exit
  - Return to Idle After Capture  Dazzle Protection
  - Return to Idle During Delay  Exposure Protection
- Default Idle Positions:**
  - Default Idle Positions
  - 550 nm
  - White
- Buttons:** Add, Copy, Remove, Test
- Time Delay:**
  - None  Manual  Automatic: 0.1 Sec.
  - Delay Position: Pre-Exposure
- Filter-Shutter / IO/LED Device:**
  - Filters:**

CoolLED-pE Filter-A	365
CoolLED-pE A	20.0 %
CoolLED-pE Filter-B	470
CoolLED-pE B	20.0 %
CoolLED-pE Filter-C	550
CoolLED-pE C	20.0 %
CoolLED-pE Filter-D	635
CoolLED-pE D	11.0 %
  - Shutters:**

CoolLED-pE Shutter-A	High
CoolLED-pE Shutter-B	High
CoolLED-pE Shutter-C	High
CoolLED-pE Shutter-D	High

- 1 Filter Group**  
Select Default Idle Positions
- 2 Copy Filter Group**  
Click Copy, enter White and click OK
- 3 Filter Group**  
Select White
- 4 Filter Controls**  
Enter the filter settings for white light from the table below
- 5 Shutter Controls**  
Set Shutters A-D to High (Open)

White Light Settings

Channel	Wavelength	Intensity %	Shutter
A	365 nm	20%	High
B	470 nm	20%	High
C	550 nm	20%	High
D	635 nm	11%	High